

**AB 5360****Rev. 1****OenoFoss™ Go****Grape Must**

This Application Brief describes the product type 'grape must' and selected parameters for this product type, available with the OenoFoss Go.

**Grape Must**

In the wine industry, the understanding of grape must is the juice from freshly crushed berries including skin, seed and pulp materials (all solids called the pomace). The must application for OenoFoss Go require a clean juice thus the pomace need to be removed from the must prior to analysis. Thus the must should be filtered.

## OenoFoss Go

OenoFoss Go is part of the FOSS wine analysis family. The main technique used for all members in the family is infrared spectroscopy – The FTIR technology. Through mathematical modelling of the infrared spectra's, the concentration of wine constituents can be determined. For each parameter in OenoFoss Go there is a prediction model calculating the results. The parameters included for grape must are: ethanol, glucose plus fructose, malic acid, total acidity and volatile acidity.

The prediction models listed in the overview below has been developed to match the OenoFoss Go software version 1.0.1.1 and newer versions.

### Prediction Model Overview

Parameter	Prediction Model							Comments
	Name	P/N	Version	Units	Performance			
					Range	Rep. <sup>1)</sup>	Acc. <sup>2)</sup>	
Total soluble solids (Brix)	Brix Must	60087603	1.0.0.1	g/100g	12-27	0.1	0.45	Total Soluble Solids are based of refraction index and calculated to common expressions - ie. Brix.
Malic acid	MalicAcid Must	60087604	1.0.0.1	g/l	1-8	0.1	0.4	
Total acidity	TotalAcidity pH 7.0 H2SO4Must	60087605	1.0.0.1	g/l	2-12	0.05	0.3	Expressed as tartaric acid
	TotalAcidity pH7.0 Tartaric Must	60089167	1.0.0.1	g/l	3-18	0.08	0.46	Expressed as tartaric acid
Yeast assimilable nitrogen	YeastAssimilableNitrogen Must	60087607	1.0.0.1	-	Low; Medium; High	NA	NA	Level indications. Medium correspond to approximative 200-300 mg/l. Low is approximate below. High is approximate beyond.
<sup>1)</sup> The absolute repeatability in corresponding unit. The repeatability expressed as pooled standard deviation of replicate. <sup>2)</sup> The absolute accuracy in corresponding unit. The accuracy is expressed as Standard Error of Prediction (SEP).								

Table 1 Overview of prediction model information

### ANALYTICS BEYOND MEASURE

**FOSS Analytical A/S**  
 Foss Allé 1  
 DK-3400 Hillerød  
 Denmark  
 Phone +45 70 10 33 70  
 Fax +45 70 10 33 71  
 E-mail info@foss.dk  
 Web www.fossanalytics.com

**FOSS Analytical Co., Ltd.**  
 6 Louyang Road, Building 1  
 215121, SIP, Suzhou  
 P.R. China  
 Phone +86 512 62 92 01 00  
 Fax +86 512 62 80 56 30  
 E-mail info@foss.dk  
 Web www.fossanalytics.com