

BENTOSTAB™ GRANULES

CLARIFICATION - FINING OF MUSTS

De-proteinizing bentonite

↓ OENOLOGICAL APPLICATIONS

BENTOSTAB™ GRANULES are activated calcium bentonite in granule form. The agglomerated bentonite offers the following advantages:

- Simpler and safer in use
- Better dispersion
- Reduced expansion times
- Increased effectiveness.

↓ INSTRUCTIONS FOR USE

Gradually add to 20 times its volume of cold water.
Leave to swell for 6 to 12 hours and add via a remontage.

↓ DOSE RATE

- 30 g to 100 g/hL on must.
- 10 g à 70 g/hL on wine.

↓ PACKAGING AND STORAGE

- 1 kg, 5 kg, 25 kg bags. **BENTOSTAB™** is marked on the packaging.

Store in a dry, odourless environment which is well ventilated at a temperature of between 5 and 25 °C.

The recommended use by date of